

BLOOD GROUP

Dr.B. Syam Sundara rao
Professor
Dept. of Pathology

BLOOD GROUP

1. who discovered the blood group system for the first time

Ans : Karl Landsteiner

Karl Landsteiner described ABO blood grouping in 1900 while Rh system was discovered later in 1940 by Landsteiner and Weiner

BLOOD GROUP

2.what is Landsteiner Law

If an agglutinogen is **present** on red cell membrane the corresponding agglutinin must **absent** in the plasma

If agglutinogen is **absent** on red cell membrane the corresponding agglutinin must **present** in the plasma

BLOOD GROUP

3. Name a few blood group systems

Ans : ABO system, Rhesus system, MNS system, Lewis system, Kelly , Duffy, Lutheran, Kidd system, Duffy system, Li and P etc.

BLOOD GROUP

4 . What type of antibodies are found in blood group “A”
?



Ans: Blood group ‘A’ has anti -B antibodies in the serum

BLOOD GROUP

5 . What type of antibodies are found in Blood group 'B'

Ans : Blood group 'B' has anti –A antibodies in the serum

ABO blood grouping system



Blood group A

If you belong to the blood group A, you have A antigens on the surface of your RBCs and B antibodies in your blood plasma.

Blood group B

If you belong to the blood group B, you have B antigens on the surface of your RBCs and A antibodies in your blood plasma.

BLOOD GROUP

6 . What type of antibodies are found in Blood group 'O'

Ans: Blood group 'O' has both Anti -A and Anti -B Antibodies in the serum

BLOOD GROUP


7. What type of antibodies are found in Blood group 'O'

Ans: Blood group 'O' has both Anti -A and Anti -B Antibodies in the serum

BLOOD GROUP


8. What type of antibodies are found in Blood group 'AB'

Ans: No antibodies



Blood group AB

If you belong to the blood group AB, you have both A and B antigens on the surface of your RBCs and no A or B antibodies at all in your blood plasma.



Blood group O

If you belong to the blood group O, you have neither A or B antigens on the surface of your RBCs but you have both A and B antibodies in your blood plasma.

Major ABO Blood Group

ABO Group	Antigen Present	Antigen Missing	Antibody Present
A	A	B	Anti-B
B	B	A	Anti-A
O	None	A and B	Anti-A&B
AB	A and B	None	None

BLOOD GROUP

9. Which type of blood group individuals are known as universal donors

Ans : Blood group 'O' individuals are known as universal donors

BLOOD GROUP

10 . Who is a universal recipient

Ans: Blood group AB individuals

BLOOD GROUP

11 . What are the methods used for blood grouping

Ans:

slide or tile method

Tube method

Microplate method

Automated method



12. What is the procedure of Tile method for blood grouping?

Procedure:

Place 1 drop of anti A and 1 drop of anti B reagent, and D reagent separately on a labelled slide or tile.

Add 1 drop of blood to each drop of the type typing antiserum.

The cells and reagent using a clean stick spread each mixture evenly on the slide or Tile over area of 10-15mm diameter.

Tilt the slide and leave the test for 2 min at room temperature (22-24) and look for agglutination.

13. What is procedure of tube method for ABO grouping?

Cell grouping:

Prepare a 2–5% cell suspension in saline from unknown blood sample.

Take 3 test tubes 1, 2, 3 and put a drop of anti-A (blue), anti-B (yellow) and anti-AB serum (pink) to them.

Add one drop of red cell suspension in each test tube.

Centrifuge at 1500rpm for one minute.

Look for agglutination either with naked eye or under the microscope.

Serum grouping:

Allow unknown blood sample stand for some time and separate the serum.

Add 2 drops of unknown serum in test tubes 4, 5, 6.

Add 1 drop of 2–5% cell suspension of known blood of A, B and O group into these test tubes.

Centrifuge at 1500rpm for one minute.

Look for agglutination either with naked eye or under the microscope.

14. What is the procedure of tube method for RH grouping?

Prepare a 2–5% cell suspension in saline from the blood to be tested.

Take a test tube and put a drop of anti-D serum

Add one drop of red cell suspension in the tube

Centrifuge at 1500rpm for 1 minute

Look for agglutination either with naked eye or under the microscope

BLOOD GROUP

15. For blood grouping which method is better tube or slide method

Ans: The tube technique is better because weaker antigens like A₂ also detected

BLOOD GROUP

16. What are the two ways for ABO blood grouping

Ans : Cell grouping (forward or direct grouping) and serum grouping (indirect or reverse grouping)

BLOOD GROUP

17. What is cell grouping or forward or direct grouping

Ans: Cell grouping or forward grouping includes the testing of antigens present on the RBC s of the patient with the help of different commercially available antisera.

BLOOD GROUP

18 . What is serum grouping or reverse or indirect grouping

Ans : Serum grouping or reverse grouping is the testing of antibodies present in the serum of the patients with the help of known cells of different blood groups

BLOOD GROUP

19 . what is cross matching

Ans: It is direct matching of blood components (serum & cell) of recipient and donor

20. What is cross matching procedure?

In a small test tube, place a drop of recipient's serum

Add washed donor red cells suspended in 5% saline

Centrifuge it at 3000rpm for one minute.

After dislodging the cell pellet gently from centrifuge tube, examine it for presence or absence of agglutination and haemolysis, first grossly and then under the low power of the microscope.

BLOOD GROUP


21 . what is the purpose of cross matching or compatibility test

Ans: The test is to prevent transfusion reactions by detecting Antibodies in recipient serum, to detect ABO incompatibilities between donor and recipient

BLOOD GROUP

22. Uses of blood grouping

Ans : For blood transfusion
In paternity disputes
Hemolytic disease of new born
Medicolegal use



BLOOD GROUP

23. Uses of blood grouping

Ans : Immunology, genetics and anthropology

Susceptibility to various diseases. Ex; o group-peptic ulcer

Part of Health checkup and job, driving licensing

BLOOD GROUP

24 . What are the adverse effects of transfusions

Ans:

Immediate

Delayed

Haemolytic(tachycardia,tachypnaea headache, hypotension ,shock)

Febrile reactions(fever, headache ,chills)

Allergic(itching, urticarial reactions, fever and chills)

BLOOD GROUP

25. What are all investigations that are performed on a donor before drawing the blood

Ans: Hb estimation, Hepatitis-B, Hepatitis c
Syphilis, HIV and malarial parasite

BLOOD GROUP

26. what do you mean by Rh incompatibility

A: Rh incompatibility occurs when Rh-ve mother has Rh + ve foetus.

The foetal blood crosses the placenta in to mothers blood to produce Ig G antibodies.

These antibodies crosses later to foetus and cause immune destruction in fetal RBC



BLOOD GROUP

27. what is bombay phenotype

Ans: The Bombay blood group was first discovered in 1952 by Bhende and others.

The phenotype is characterised by the absence of A, B, H antigens on the red cells. The serum of these persons contain anti A, anti B, and Anti H antibodies


BLOOD GROUP

28. what is coombs(Anti globulin) test

Ans: It is used to demonstrate presence of Ig G antibodies against red blood cells in the patient
IgG antibodies do not cause agglutination of RBCs, coombs serum is added which combines with IgG on RBCs surface resulting agglutination

BLOOD GROUP

29.what are Direct combs test applications

- ▶ A:Hemolytic transfusion reactions
 - ▶ Auto hemolytic anemia
 - ▶ Hemolytic disease of new born
 - ▶ Drug induced hemolysis
- 

BLOOD GROUP

30.what are Indirect coombs test applications

Indirect:

- ▶ Compatibility testing
- ▶ Detection of unexpected antibody(IgG)s to Rh factor in serum of pregnant patients
- ▶ Detection of red cell antigens not detected by other technique

BLOOD GROUP

31.what are blood components

Ans: A Single unit of whole blood can be separated in to one unit of packed red cells, one unit of platelets and one unit of fresh frozen plasma

32. What is one unit of blood?

- ▶ **ANS:** 350ml of blood in 49ml of anticoagulant citrate phosphate–Dextrose Adenine

33. What is the difference between complete and incomplete antibodies?

- ▶ ANS: Naturally occurring anti A, anti B antibodies are IGM types– complete antibodies and cause agglutination in saline suspension.
- ▶ Incomplete antibodies IgG type and result because of antigenic stimulation and don't cause agglutination in saline suspension.

34. What test is employed to detect low level antibodies against red cells?

- ▶ ANS: Coomb tests.

35. Which are the immune haemolytic anemias in which coombs test is indicated?

- ▶ ANS: Auto immune haemolytic anemia SLE, cold agglutinin disease, paroxysmal nocturnal haemoglobinuria, infectious mononucleosis.

36. How are blood and its components stored?

- ▶ ANS: Blood is stored in a refrigerator at 2–4 c upto 4–6 weeks, and platelet can be stored up to 5 days .

THANK YOU

